

- Single use
- Cost-effective
- Reduces cross infection
- Latex free
- Colour coded sizes



Single Use Tourniquet Cuffs

Denyer Disposable Tourniquet Cuffs are sterile, single-use and come in a variety of sizes, colour coded for easy identification.

All cuffs are available with single or dual hose with conventional snaplock fitting and most are also available with Leur Lock connection.

As well as our disposable cuffs, we carry a comprehensive range of reusable cuffs, tourniquet machines and accessories.

Tourniquet Cuff Sleeve

Denyers Limb protection sleeves are specifically designed to protect patients limbs from possible wrinkling, pinching and shearing of skin and soft tissues of the limb.

Each Limb Protection Sleeve is intended for use with a specific Low pressure Tourniquet Cuff.

Within the range of limb sizes recommended for the matching Cuff, the sleeve helps to protect while applying a low pressure that is less than a snugly wrapped Cuff.



Specifications - Disposable Cuffs

Cuff Sleeves



Single Disposable Cuff



Item Code	Size	Description	Unit	Sleeve Code	Unit
TCQ/D082-1	Infant 8"	Single hose Light Blue	Box 10	TCQ/9002	Box of 40
TCQ/D123-1	Child 12"	Single hose Burgundy	Box 10	TCQ/9003	Box of 40
TCQ/D184-1	Std Arm 18"	Single hose Red	Box 10	TCQ/9003	Box of 40
TCQ/D184-1L	Std Arm 1	8" Single hose Red Luer	Box 10	TCQ/9003	Box of 40
TCQ/D244-1	Std Leg 24"	Single hose Yellow	Box 10	TCQ/9004	Box of 40
TCQ/D344-1	Long Leg 34"	Single hose Blue	Box 10	TCQ/9004	Box of 40
TCQ/D344-1L	Long Leg 34"	Single hose Blue Luer	Box 10	TCQ/9004	Box of 40
TCQ/D444-1	XLong Leg 44"	Single hose Green	Box 10	TCQ/9005	Box of 25
TCQ/D444-1L	XLong Leg 44"	Single hose Green Luer	Box 10	TCQ/9005	Box of 25
TCQ/DC344-1	Conical Long Leg 34"	Single hose	Box 10	TCQ/9004	Box of 40
TCQ/DC344-1L	Conical Long Leg 34"	Single hose Luer	Box 10	TCQ/9004	Box of 40



Double Disposable Cuff



TCQ/D082-2	Infant 8"	Double hose Light Blue	Box 10	TCQ/9002	Box of 40
TCQ/D123-2	Child 12"	Double hose Burgundy	Box 10	TCQ/9003	Box of 40
TCQ/D184-2	Std Arm 18"	Double hose Red	Box 10	TCQ/9003	Box of 40
TCQ/D244-2	Std Leg 24"	Double hose Yellow	Box 10	TCQ/9004	Box of 40
TCQ/D344-2	Long Leg 34"	Double hose Blue	Box 10	TCQ/9004	Box of 40
TCQ/D444-2	XLong Leg 44"	Double hose Green	Box 10	TCQ/9005	Box of 25

DENYERS INTERNATIONAL LTD

Freecall: 1800 DENYER (1800 336 937)

42- 46 Nelson Street, Moorabbin,
Victoria 3189 Australia
Email: sales@denyers.com.au
www.denyers.com.au

DENYERS NZ LTD

Freecall: 0800 DENYERS (0800 336 9377)

Phone: +64 9 215 6671
Fax: +64 9 9732474
Email: sales@denyers.co.nz
www.denyers.co.nz

DENYERS
INTERNATIONAL

PROUDLY AUSTRALIAN