



HERITAGE MEDICAL PRODUCTS

THE WEDGE IS OUR EDGE

Experience the superior quality of the Heritage Medical Sonobed®.



HMP is a Sonographer owned company specializing in ergonomic exam beds.

Our custom designs offer patented features that enable Sonographers to place their patient in optimal positions and obtain the best views possible. Easing strain and stress off of the Sonographer, giving the patient a comfy spot to rest during the exam.



QUALITY ASSURED

HERITAGE MEDICAL PRODUCTS

CONTENTS

HMP History.....	Page 1
Sonobed Multi-Specialty.....	Page 2
Sonobed Cardiac.....	Page 3
2000 Series Vascular.....	Page 4
9000 Series Cardiac.....	Page 5

CONTACT

At Heritage Medical Products, we are committed to answering any questions you have 24/7/365. Call any time at 1-888-888-8888 and one of our friendly technicians will assist you with questions about any Sonobed product.

Or visit our website at www.sonobed.com



2

HMP SONOBED™ MULTI-SPECIALTY FDA REGISTERED



Double Wedge Sections™
STANDARD FEATURE

Patented Double Power Wedge Sections™ are optional.
Powered for patient positioning assistance.

Ergo Side Rails
OPTIONAL FEATURE

Straps
STANDARD FEATURE

Hidden and retractable so they never leave the bed.

Central Locking Casters
OPTIONAL FEATURE

FEATURES AND BENEFITS

Double Wedge Section™ (2 - 10"x13" patient positioning Wedge Sections™)	Infinite Position Electric Hand Control (Height & Fowler)
UL Listed 9601-1 in accordance with CAN/CSA No 601	Heavy-duty Steel Frame Painted with Baked on Powder Coated Finish for Years of Durability
Electric Positioning Fowler	600 lb Load and Lift Capacity
Electric Height Adjustable from 26" to 38"	Spacious Measurements - 30" Width x 71" Length
5" Total Lock Casters	Convenient Storage Shelf
Four Protective Rubber Corner Bumpers	Custom Designed Foot Panel

WWW.SONOBED.COM




QUALITY ASSURED

HERITAGE MEDICAL PRODUCTS

*HMP has an
exclusive patent on
the Power Wedge Section™
and Power Cardiac
Drop Section combined
technology!*

3

HMP SONOBED™ CARDIAC FDA REGISTERED



Wedge Section™
STANDARD FEATURE

A BUILT-IN patient positioning wedge for
patients to lay against during the scan.

Headrail
OPTIONAL FEATURE

Power Drop Section™
OPTIONAL FEATURE

Designed for multi-hand scanning and as Scansophers
Early user and close the Drop Section with standard
hand control and/or optional foot control.
Allows for unobstructed access.

Power Foot Panel
OPTIONAL FEATURE




FEATURES AND BENEFITS

13" x 13" Imaging Drop Section - Manual	600 lb Load Capacity
13" x 13" Wedge Section™ - Manually rises up & holds your patient comfortable on their left side	Infinite Position Electric Hand Control (Height & Fowler)
Electric Positioning Fowler	Spacious Measurements - 30" Width x 71" Length
Powerful & Fast Electric Height Adjustable from 26" to 38" with a 600 lb Lift Capacity	Heavy duty Steel Frame Painted with Baked on Powder Coated Finish for Years of Durability
UL Listed 9601-1 in accordance with CAN/CSA No 601	Convenient Storage Shelf
5" Total Lock Casters with Corner Bumpers	Custom Designed Foot Panel

WWW.SONOBED.COM

4

2000 SERIES VASCULAR FDA REGISTERED



Double Wedge Sections™ STANDARD FEATURE

Patented Double Power Wedge Sections™ are control powered for patient positioning assistance.

Armboard
STANDARD

Footboard
STANDARD

Central Locking Casters
OPTIONAL FEATURE



Electronic Trendelenburg 15°/30° Reverse
STANDARD ON VASCULAR BED ONLY

FEATURES AND BENEFITS

Double Wedge Section™ (2 - 13"V3" patient positioning Wedge Sections™)	Infinite Position Hand Control (Height, Fowler & Optional Trendelenburg)
Electric Trendelenburg +/- 30° Vascular Bed only (Standard feature on our vascular models Sonobed™ and 2000 Series)	Heavy-duty Steel Frame Painted with Baked on Powder Coated Finish for Years of Durability
Electric Positioning Fowler	600 lb Load and Lift Capacity
Electric Height: Adjustable from 36" to 38"	Spacious Measurements - 30" Width x 74" Length
UL Listed 9501-1 in accordance with CAN/CSA No 601	Four Protective Rubber Corner Bumpers
5" Total Lock Casters	Convenient Storage Shelf

WWW.SONOBED.COM

5

2000 SERIES CARDIAC FDA REGISTERED



Wedge Section™ STANDARD FEATURE

A BUILT-IN patient positioning wedge for patients to lay against during the scan.

Power Drop Section™ OPTIONAL FEATURE

Designed for right hand scanning cardiac Sonographs. Easily open and close the Drop Section with standard hand control and/or optional foot control. Allows for unobstructed apical views.

12" Low Setting



49" Height Setting



Ultra Low Column OPTIONAL FEATURE

Adjustable bed height allows you to lower down to 12" or raise up to 49".

QUALITY ASSURED

HERITAGE MEDICAL PRODUCTS

We stand behind
the quality of every
bed we manufacture
to your complete
satisfaction...
GUARANTEED!

FEATURES AND BENEFITS

13" x 13" Imaging Drop Section - Manual	600 lb Load Capacity
13" x 13" Wedge Section™ - Manually rises up & holds your patient comfortable on their left side	Infinite Position Electric Hand Control (Height & Fowler)
Electric Positioning Fowler	Spacious Measurements - 30" Width x 74" Length
Powerful & Fast Electric Height Adjustable from 26" to 38" with a 600 lb Lift Capacity	Heavy-duty Steel Frame Painted with Baked on Powder Coated Finish for Years of Durability
UL Listed 9501-1 in accordance with CAN/CSA No 601	Four Protective Rubber Corner Bumpers
5" Total Lock Casters	Convenient Storage Shelf for Supplies

WWW.SONOBED.COM

4

2000 SERIES VASCULAR

FDA REGISTERED



Double Wedge Section™

STANDARD FEATURE

Patented Double Power Wedge Sections™ are optional. Powered for patient positioning assistance.



Armboard
STANDARD

Footboard
STANDARD

Central Locking Casters
OPTIONAL FEATURE

Electronic Trendelenburg 15°/30° Reverse
STANDARD ON VASCULAR BED ONLY

FEATURES AND BENEFITS

Double Wedge Section™ (2 - 19½" patient positioning Wedge Sections™)	Infinite Position Hand Control (Height, Fowler & Optional Trendelenburg)
Electric Trendelenburg +/- 30° Vascular Bed only (Standard feature on our vascular models Sonosed™ and 2000 Series)	Heavy-duty Steel Frame Painted with Baked on Powder Coated Finish for Years of Durability
Electric Positioning Fowler	600 lb Load and Lift Capacity
Electric Height: Adjustable from 26" to 38"	Spacious Measurements - 30" Width x 74" Length
UL Listed 2601-1 in accordance with CAN/CSA No. 601	Four Protective Rubber Corner Bumpers
5" Total Lock Casters	Convenient Storage Shelf

WWW.SONOBED.COM

5

2000 SERIES CARDIAC

FDA REGISTERED

Wedge Section™

STANDARD FEATURE

A BUILT IN patient positioning wedge for patients to lay against during the exam.

Power Drop Section™

OPTIONAL FEATURE

Designed for right hand identity cardiac Stereographs. Easily open and close the Drop Section with standard hand control and/or optional foot control. Allows for undistorted axial views.

12" Low Setting



49" Height Setting

Ultra Low Column
OPTIONAL FEATURE

Adjustable bed height allows you to lower down to 19" or rise up to 49".

QUALITY ASSURED

HERITAGE MEDICAL PRODUCTS

We stand behind
the quality of every
bed we manufacture
to your complete
satisfaction...
GUARANTEED!

FEATURES AND BENEFITS

13" x 13" Imaging Drop Section - Manual	600 lb Load Capacity
13" x 13" Wedge Section™ - Manually rises up & holds your patient comfortable on their left side	Infinite Position Electric Hand Control (Height & Fowler)
Electric Positioning Fowler	Spacious Measurements - 30" Width x 74" Length
Powerful & Fast Electric Height Adjustable from 26" to 38" with a 600 lb Lift Capacity	Heavy-duty Steel Frame Painted with Baked on Powder Coated Finish for Years of Durability
UL Listed 2601-1 in accordance with CAN/CSA No. 601	Four Protective Rubber Corner Bumpers
5" Total Lock Casters	Convenient Storage Shelf for Supplies

WWW.SONOBED.COM



QUALITY ASSURED

HERITAGE MEDICAL PRODUCTS



QUALITY ASSURED

HERITAGE MEDICAL PRODUCTS